## VLASTNOSTI VANKÚŠA JAY BALANCE

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## Referencie

1. Brienza, D., Kelsey, S., Karg, P., Allegretti, A., Olson, M., Schmeler, M., Holm, M. (2010). A Randomized Clinical Trial on Preventing Pressure Ulcers with Wheelchair Seat Cushions. *Journal of the American Geriatrics Society*, 58(12), 2308–2314. doi: 10.1111/j.1532-5415.2010.03168.x
2. Call, E., Hetzel, T., Mclean, C., Burton, J. N., & Oberg, C. (2017). Off loading wheelchair cushion provides best case reduction in tissue deformation as indicated by MRI. *Journal of Tissue Viability*, 26(3), 172–179. doi: 10.1016/j.jtv.2017.05.002
3. Garcia-Mendez, Y., Pearlman, J. L., Boninger, M. L., & Cooper, R. A. (2013). Health risks of vibration exposure to wheelchair users in the community. *The Journal of Spinal Cord Medicine*, 36(4), 365–375. doi: 10.1179/2045772313y.0000000124
4. Shabshin, N., Zoizner, G., Herman, A., Ougortsin, V., & Gefen, A. (2010). Use of weight-bearing MRI for evaluating wheelchair cushions based on internal soft-tissue deformations under ischial tuberosities. *The Journal of Rehabilitation Research and Development*, 47(1), 31. doi: 10.1682/jrrd.2009.07.0105
5. Watanabe, L. (2017, September). Immersion, envelopment and off-loading. *Mobility Management*. Retrieved from <https://mobilitymgmt.com/Articles/2017/09/01/Seating-Strategies.aspx>

DISCLAIMER: FOR PROFESSIONAL USE ONLY. THIS DOCUMENT (AND THE WEBSITE REFERENCED HEREIN) DO NOT PROVIDE MEDICAL ADVICE. Sunrise Medical does not provide clinician services. The information contained on this document (and the website referenced herein), including, but not limited to, the text, graphics, images, and descriptions, are for informational purposes only and should be utilized as a general resource for clinicians and suppliers to then use clinical reasoning skills to determine optimal seating and mobility solutions for individual patients. No material on this document (or on the website) is intended to be used as (or a substitute for) professional medical advice, diagnosis or treatment. Clinicians should adhere to their professional medical training and their institution's or practice's recommended guidelines. Reliance on this document (and the information contained herein) is solely at your own risk.